

Zymo-Teck® the secret of quality grafts



Bioteck bone substitutes are obtained from equine bone tissue treated with Zymo-Teck®. This exclusive proprietary process is based on the utilization of lytic enzymes operating at controlled temperatures. This enables the complete elimination of the antigen components of the tissue, without the mineral phase undergoing any changes.



The unmodified bone mineral component is recognized as endogenous by the osteoclasts, thereby allowing for the total remodeling of the graft, which is completely replaced, in physiological time, by new patient vital bone tissue. The best possible condition for osseointegrated implants.



enzymatic process



beta ray sterilization



safety and quality



total biocompatibility




complete remodelling



Bioteck S.p.A.
 Headquarters:
 Via E. Fermi 49 - 36057 Arcugnano (Vicenza) - Italy
 Tel. +39 0444 289366 - fax: +39 0444 285272
 info@bioteck.com - www.bioteck.com
 Production and R&D Center:
 Via G. Agnelli, 3 - 10020 Riva presso Chieri (Turin) - Italy

bioteck.com



OSTEOPLANT®
BIOCOLLAGEN®
BIO-GEN®

Natural Bone Substitutes
 & Membranes



BIO-GEN®

Granules

| | | | | |
|------------------|----------------------------------|-------|-------------|----------|
| BGS-05s* | Cancellous Granules | 1 btl | 0.5g=1cc | 0.25-1mm |
| BGS-09s* | Cancellous Granules | 1 btl | 0.5g=1cc | 1-2mm |
| BGS-11s* | Cancellous Granules | 1 btl | 1g=2cc | 1-2mm |
| BGS-23s* | Cancellous Granules | 1 btl | 1g=2cc | 2-3mm |
| BGS-20 | Cancellous Granules | 1 btl | 2g=0,25 | 1mm |
| BGS-22 | Cancellous Granules | 1 btl | 2g=4cc | 1-2mm |
| BGC-05s* | Cortical Granules | 1 btl | 0.5g=1cc | 0.25-1mm |
| BGM-10s* | Cancellous Cortical Mix Granules | 1 btl | 0.25g=0.5cc | 0.25-1mm |
| BGM-05s* | Cancellous Cortical Mix Granules | 1 btl | 0.5g=1cc | 0.25-1mm |
| BGM-100s* | Cancellous Cortical Mix Granules | 1 btl | 1g=2cc | 0.25-1mm |
| BGM-20 | Cancellous Cortical Mix Granules | 1 btl | 2g=4cc | 0.5-1mm |

*: The product is in single sterile packaging



BIO-GEN®

Putty

| | | | |
|-----------------|----------------------|------|-------|
| BGP-01s* | Cancellous Dry Paste | 1 pc | 0.5cc |
|-----------------|----------------------|------|-------|

*: The product is in single sterile packaging



BIO-GEN®

Cancellous Blocks and Wedge

| | | | |
|---------------|-----------------------|------|-----------------|
| BGB-11 | Cancellous Bone Block | 1 pc | 10 x 10 x 10 mm |
| BGB-12 | Cancellous Bone Block | 1 pc | 10 x 10 x 20 mm |
| BGB-30 | Cancellous Bone Wedge | 1 pc | 25 x 10 x 5 mm |



BIOCOLLAGEN®

Collagen Membrane

| | | | |
|-----------------|-------------------|------|------------------|
| BCG-01s* | Collagen Membrane | 1 pc | 25 x 25 x 0.2 mm |
| BCG-02s* | Collagen Membrane | 1 pc | 15 x 20 x 0.2 mm |

*: The product is in single sterile packaging

GO TO DYNAMIC LITERATURE


Bio-Gen



Biocollagen



Osteoplant



BIO-GEN®

GEL Granules

| | | | |
|--------------------|-----------------------------|----------------|---------|
| BGM-GEL05s* | Cancellous Cortical Mix Gel | 1 syr / 0.5ml | 0.5-1mm |
| BGM-GEL1s* | Cancellous Cortical Mix Gel | 1 syr / 1ml | 0.5-1mm |
| BGM-GEL02n* | Cancellous Cortical Mix Gel | 1 syr / 0,25ml | 0,5-1mm |

*: The product is in single sterile packaging



OSTEOPLANT®

Demineralized Bone Matrix (DBM)

| | | | |
|--------------------|------------------------------|-------|--------|
| OGS-AC5n* | Activagen DBM granules | 1 btl | 0.5 cc |
| OGS-ACM501* | Activagen DBM moldable paste | 1 syr | 0.5 ml |
| OGS-ACM601* | Activagen DBM moldable paste | 1 syr | 1 ml |
| OGS-GEL1n* | Angiostad DBM gel | 1 syr | 1 ml |

*: The product is in single sterile packaging